

STRONG ACID AND ALKALI STORAGE CABINETS



STRONG ACID AND ALKALI STORAGE CABINETS

The Strong Acid and Alkali Cabinet is robust, purpose-built storage unit which is potentially designed to store strong and reactive acids and alkalis. This is a containment of highly corrosive substances with the capacity that ranges from 45L- 340L offers ample space and has separate shelves for organised storage. It reduces the risk of accidental spills, exposure, or reactions and ensuring the safety of personnel and the surrounding environment.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

LTCAB9 STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

Model	LTCAB9-1	LTCAB9-2	LTCAB9-3	LTCAB9-4
Capacity	4 gallon (15L)	12 gallon (45L)	20 gallon (75L)	28 gallon (105L)
Material	pp			
Adjustable Shelf	Optional	1 pc		2 pcs
Door	Single / Manual		Double / Manual	
External Size	430x430x500 mm	590x460x890 mm	910x590x890 mm	910x600x900 mm
Package Size	515x515x685 mm	675x545x1025 mm	995x675x1035 mm	995x685x1035 mm
Gross Weight	16 Kg	30 Kg	46 Kg	48 Kg


Model	LTCAB9-5	LTCAB9-6	LTCAB9-7	LTCAB9-8
Capacity	30 gallon (114L)	45 gallon (170L)	60 gallon (227L)	90 gallon (340L)
Material	pp			
Adjustable Shelf	2 pcs			
Door	Double / Manual			
External Size	1090x460x1120mm	1090x460x1650mm	860x860x1650 mm	1090x860x1650 mm
Package Size	1175x545x1255 mm	1175x545x1785 mm	945x945x1785 mm	1173x943x1785 mm
Gross Weight	58 Kg	72 Kg	90 Kg	107 Kg




LTCAB9-1




LTCAB9-2



LTCAB9-3



LTCAB9-4



LTCAB9-5



LTCAB9-6



LTCAB9-7



LTCAB9-8

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com